

# Bluetooth<sup>®</sup> Low Energy Protocol Stack

API Reference Manual: BLP

Renesas MCU

Target Device

RL78/G1D

All information contained in these materials, including products and product specifications, represents information on the product at the time of publication and is subject to change by Renesas Electronics Corp. without notice. Please review the latest information published by Renesas Electronics Corp. through various means, including the Renesas Electronics Corp. website (<http://www.renesas.com>).

## Notice

1. Descriptions of circuits, software and other related information in this document are provided only to illustrate the operation of semiconductor products and application examples. You are fully responsible for the incorporation of these circuits, software, and information in the design of your equipment. Renesas Electronics assumes no responsibility for any losses incurred by you or third parties arising from the use of these circuits, software, or information.
2. Renesas Electronics has used reasonable care in preparing the information included in this document, but Renesas Electronics does not warrant that such information is error free. Renesas Electronics assumes no liability whatsoever for any damages incurred by you resulting from errors in or omissions from the information included herein.
3. Renesas Electronics does not assume any liability for infringement of patents, copyrights, or other intellectual property rights of third parties by or arising from the use of Renesas Electronics products or technical information described in this document. No license, express, implied or otherwise, is granted hereby under any patents, copyrights or other intellectual property rights of Renesas Electronics or others.
4. You should not alter, modify, copy, or otherwise misappropriate any Renesas Electronics product, whether in whole or in part. Renesas Electronics assumes no responsibility for any losses incurred by you or third parties arising from such alteration, modification, copy or otherwise misappropriation of Renesas Electronics product.
5. Renesas Electronics products are classified according to the following two quality grades: "Standard" and "High Quality". The recommended applications for each Renesas Electronics product depends on the product's quality grade, as indicated below.  
"Standard": Computers; office equipment; communications equipment; test and measurement equipment; audio and visual equipment; home electronic appliances; machine tools; personal electronic equipment; and industrial robots etc.  
"High Quality": Transportation equipment (automobiles, trains, ships, etc.); traffic control systems; anti-disaster systems; anti-crime systems; and safety equipment etc.  
  
Renesas Electronics products are neither intended nor authorized for use in products or systems that may pose a direct threat to human life or bodily injury (artificial life support devices or systems, surgical implantations etc.), or may cause serious property damages (nuclear reactor control systems, military equipment etc.). You must check the quality grade of each Renesas Electronics product before using it in a particular application. You may not use any Renesas Electronics product for any application for which it is not intended. Renesas Electronics shall not be in any way liable for any damages or losses incurred by you or third parties arising from the use of any Renesas Electronics product for which the product is not intended by Renesas Electronics.
6. You should use the Renesas Electronics products described in this document within the range specified by Renesas Electronics, especially with respect to the maximum rating, operating supply voltage range, movement power voltage range, heat radiation characteristics, installation and other product characteristics. Renesas Electronics shall have no liability for malfunctions or damages arising out of the use of Renesas Electronics products beyond such specified ranges.
7. Although Renesas Electronics endeavors to improve the quality and reliability of its products, semiconductor products have specific characteristics such as the occurrence of failure at a certain rate and malfunctions under certain use conditions. Further, Renesas Electronics products are not subject to radiation resistance design. Please be sure to implement safety measures to guard them against the possibility of physical injury, and injury or damage caused by fire in the event of the failure of a Renesas Electronics product, such as safety design for hardware and software including but not limited to redundancy, fire control and malfunction prevention, appropriate treatment for aging degradation or any other appropriate measures. Because the evaluation of microcomputer software alone is very difficult, please evaluate the safety of the final products or systems manufactured by you.
8. Please contact a Renesas Electronics sales office for details as to environmental matters such as the environmental compatibility of each Renesas Electronics product. Please use Renesas Electronics products in compliance with all applicable laws and regulations that regulate the inclusion or use of controlled substances, including without limitation, the EU RoHS Directive. Renesas Electronics assumes no liability for damages or losses occurring as a result of your noncompliance with applicable laws and regulations.
9. Renesas Electronics products and technology may not be used for or incorporated into any products or systems whose manufacture, use, or sale is prohibited under any applicable domestic or foreign laws or regulations. You should not use Renesas Electronics products or technology described in this document for any purpose relating to military applications or use by the military, including but not limited to the development of weapons of mass destruction. When exporting the Renesas Electronics products or technology described in this document, you should comply with the applicable export control laws and regulations and follow the procedures required by such laws and regulations.
10. It is the responsibility of the buyer or distributor of Renesas Electronics products, who distributes, disposes of, or otherwise places the product with a third party, to notify such third party in advance of the contents and conditions set forth in this document, Renesas Electronics assumes no responsibility for any losses incurred by you or third parties as a result of unauthorized use of Renesas Electronics products.
11. This document may not be reproduced or duplicated in any form, in whole or in part, without prior written consent of Renesas Electronics.
12. Please contact a Renesas Electronics sales office if you have any questions regarding the information contained in this document or Renesas Electronics products, or if you have any other inquiries.

(Note 1) "Renesas Electronics" as used in this document means Renesas Electronics Corporation and also includes its majority-owned subsidiaries.

(Note 2) "Renesas Electronics product(s)" means any product developed or manufactured by or for Renesas Electronics.

## General Precautions in the Handling of MPU/MCU Products

The following usage notes are applicable to all MPU/MCU products from Renesas. For detailed usage notes on the products covered by this document, refer to the relevant sections of the document as well as any technical updates that have been issued for the products.

### 1. Handling of Unused Pins

Handle unused pins in accordance with the directions given under Handling of Unused Pins in the manual.

- The input pins of CMOS products are generally in the high-impedance state. In operation with an unused pin in the open-circuit state, extra electromagnetic noise is induced in the vicinity of LSI, an associated shoot-through current flows internally, and malfunctions occur due to the false recognition of the pin state as an input signal become possible. Unused pins should be handled as described under Handling of Unused Pins in the manual.

### 2. Processing at Power-on

The state of the product is undefined at the moment when power is supplied.

- The states of internal circuits in the LSI are indeterminate and the states of register settings and pins are undefined at the moment when power is supplied.  
In a finished product where the reset signal is applied to the external reset pin, the states of pins are not guaranteed from the moment when power is supplied until the reset process is completed. In a similar way, the states of pins in a product that is reset by an on-chip power-on reset function are not guaranteed from the moment when power is supplied until the power reaches the level at which resetting has been specified.

### 3. Prohibition of Access to Reserved Addresses

Access to reserved addresses is prohibited.

- The reserved addresses are provided for the possible future expansion of functions. Do not access these addresses; the correct operation of LSI is not guaranteed if they are accessed.

### 4. Clock Signals

After applying a reset, only release the reset line after the operating clock signal has become stable. When switching the clock signal during program execution, wait until the target clock signal has stabilized.

- When the clock signal is generated with an external resonator (or from an external oscillator) during a reset, ensure that the reset line is only released after full stabilization of the clock signal. Moreover, when switching to a clock signal produced with an external resonator (or by an external oscillator) while program execution is in progress, wait until the target clock signal is stable.

### 5. Differences between Products

Before changing from one product to another, i.e. to a product with a different part number, confirm that the change will not lead to problems.

- The characteristics of an MPU or MCU in the same group but having a different part number may differ in terms of the internal memory capacity, layout pattern, and other factors, which can affect the ranges of electrical characteristics, such as characteristic values, operating margins, immunity to noise, and amount of radiated noise. When changing to a product with a different part number, implement a system-evaluation test for the given product.

# How to Use This Manual

## 1. Purpose and Target Readers

This manual describes the API (Application Program Interface) of the Blood Pressure profile (BLP) of the Bluetooth Low Energy protocol stack (BLE software), which is used to develop Bluetooth applications that incorporate the Renesas Bluetooth low energy microcontroller RL78/G1D. It is intended for users designing application systems incorporating this software. A basic knowledge of microcontrollers and Bluetooth low energy is necessary in order to use this manual.

### Related documents

The related documents indicated in this publication may include preliminary versions. However, preliminary versions are not marked as such.

Document Name	Document No.
Bluetooth Low Energy Protocol Stack	
User's Manual	R01UW0095E
API Reference Manual: Basics	R01UW0088E
API Reference Manual: FMP	R01UW0089E
API Reference Manual: PXP	R01UW0090E
API Reference Manual: HTP	R01UW0091E
API Reference Manual: BLP	This manual
API Reference Manual: HOGP	R01UW0093E
API Reference Manual: ScPP	R01UW0094E
API Reference Manual: HRP	R01UW0097E
API Reference Manual: CSCP	R01UW0098E
API Reference Manual: CPP	R01UW0099E
API Reference Manual: GLP	R01UW0103E
API Reference Manual: TIP	R01UW0106E
API Reference Manual: RSCP	R01UW0107E
API Reference Manual: ANP	R01UW0108E
API Reference Manual: PASP	R01UW0109E
API Reference Manual: LNP	R01UW0113E
Application Note: Sample Program	R01AN1375E
Application Note: rBLE Command Specification	R01AN1376E

## List of Abbreviations and Acronyms

Abbreviation	Full Form	Remark
ANP	Alert Notification Profile	
ANS	Alert Notification Service	
API	Application Programming Interface	
ATT	Attribute Protocol	
BAS	Battery Service	
BB	Base Band	
BD_ADDR	Bluetooth Device Address	
BLE	Bluetooth low energy	
BLP	Blood Pressure Profile	
BLS	Blood Pressure Service	
CPP	Cycling Power Profile	
CPS	Cycling Power Service	
CSCP	Cycling Speed and Cadence Profile	
CSCS	Cycling Speed and Cadence Service	
CSRK	Connection Signature Resolving Key	
CTS	Current Time Service	
DIS	Device Information Service	
EDIV	Encrypted Diversifier	
FMP	Find Me Profile	
GAP	Generic Access Profile	
GATT	Generic Attribute Profile	
GLP	Glucose Profile	
GLS	Glucose Service	
HCI	Host Controller Interface	
HID	Human Interface Device	
HIDS	HID Service	
HOGP	HID over GATT Profile	
HRP	Heart Rate Profile	
HRS	Heart Rate Service	
HTP	Health Thermometer Profile	
HTS	Health Thermometer Service	
IAS	Immediate Alert Service	
IRK	Identity Resolving Key	
L2CAP	Logical Link Control and Adaptation Protocol	
LE	Low Energy	

Abbreviation	Full Form	Remark
LL	Link Layer	
LLS	Link Loss Service	
LNP	Location and Navigation Profile	
LNS	Location and Navigation Service	
LTK	Long Term Key	
MCU	Micro Controller Unit	
MITM	Man-in-the-middle	
MTU	Maximum Transmission Unit	
NDCS	Next DST Change Service	
OOB	Out of Band	
OS	Operating System	
PASP	Phone Alert Status Profile	
PASS	Phone Alert Status Service	
PXP	Proximity Profile	
RF	Radio Frequency	
RSCP	Running Speed and Cadence Profile	
RSCS	Running Speed and Cadence Service	
RSSI	Received Signal Strength Indication	
RTUS	Reference Time Update Service	
ScPP	Scan Parameters Profile	
ScPS	Scan Parameters Service	
SM	Security Manager	
SMP	Security Manager Protocol	
STK	Short Term Key	
TIP	Time Profile	
TK	Temporary Key	
TPS	Tx Power Service	
UART	Universal Asynchronous Receiver Transmitter	
UUID	Universal Unique Identifier	

Abbreviation	Full Form	Remark
APP	Application	
CSI	Clocked Serial Interface	
IIC	Inter-Integrated Circuit	
RSCIP	Renesas Serial Communication Interface Protocol	
VS	Vendor Specific	

All trademarks and registered trademarks are the property of their respective owners.

Bluetooth is a registered trademark of Bluetooth SIG, Inc. U.S.A.

EEPROM is a trademark of Renesas Electronics Corporation.

Windows, Windows NT and Windows XP are registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

PC/AT is a trademark of International Business Machines Corporation.

# Contents

1. Overview.....	1
2. Common Definitions .....	2
2.2 Status Definitions.....	4
3. Blood Pressure Profile.....	5
3.1 Definitions .....	5
3.2 Functions .....	12
3.2.1 RBLE_BLP_Sensor_Enable.....	13
3.2.2 RBLE_BLP_Sensor_Disable.....	13
3.2.3 RBLE_BLP_Sensor_Send_Measurements.....	14
3.2.4 RBLE_BLP_Collector_Enable.....	15
3.2.5 RBLE_BLP_Collector_Disable.....	16
3.2.6 RBLE_BLP_Collector_Read_Char .....	17
3.2.7 RBLE_BLP_Collector_Write_Char .....	17
3.3 Events .....	18
3.3.1 RBLE_BLP_EVENT_SENSOR_ENABLE_COMP.....	19
3.3.2 RBLE_BLP_EVENT_SENSOR_DISABLE_COMP.....	19
3.3.3 RBLE_BLP_EVENT_SENSOR_ERROR_IND .....	19
3.3.4 RBLE_BLP_EVENT_SENSOR_SEND_MEASUREMENTS_COMP.....	19
3.3.5 RBLE_BLP_EVENT_SENSOR_CFG_INDNTF_IND .....	20
3.3.6 RBLE_BLP_EVENT_SENSOR_COMMAND_DISALLOWED_IND .....	20
3.3.7 RBLE_BLP_EVENT_COLLECTOR_ENABLE_COMP.....	21
3.3.8 RBLE_BLP_EVENT_COLLECTOR_DISABLE_COMP .....	22
3.3.9 RBLE_BLP_EVENT_COLLECTOR_ERROR_IND .....	22
3.3.10 RBLE_BLP_EVENT_COLLECTOR_MEASUREMENTS_IND .....	23
3.3.11 RBLE_BLP_EVENT_COLLECTOR_READ_CHAR_RESPONSE.....	23
3.3.12 RBLE_BLP_EVENT_COLLECTOR_WRITE_CHAR_RESPONSE.....	24
3.3.13 RBLE_BLP_EVENT_COLLECTOR_COMMAND_DISALLOWED_IND .....	24
3.4 Message Sequence Chart .....	25



4. Notes .....	26
Appendix A How to Read Definition Tables.....	27
Appendix B Referenced Documents .....	29
Appendix C Terminology .....	30

## 1. Overview

This manual describes the API (Application Program Interface) of the Blood Pressure profile (BLP) of the Bluetooth Low Energy protocol stack (BLE software), which is used to develop Bluetooth applications that incorporate Renesas Bluetooth low energy microcontroller RL78/G1D.

For details about the organization and features of BLE software, see the Bluetooth Low Energy Protocol Stack User's Manual.

## 2. Common Definitions

This section describes the definitions common to the API of each profile.

### 2.1 Service Definitions

This section describes the common definitions of services used by the API of multiple profiles.

Declaration of enumerated type for alert level

```
enum RBLE_SVC_ALT_LVL_enum {
    RBLE_SVC_ALERT_NONE = 0x00,           No alert
    RBLE_SVC_ALERT_MILD,                   Mild alert
    RBLE_SVC_ALERT_HIGH                    High alert
};
```

Declaration of enumerated type for PnP ID characteristic vendor ID field

```
enum RBLE_SVC_PNP_VENDOR_ID_enum {
    RBLE_SVC_SIG_ASSIGNED_ID = 0x01,      Vendor ID assigned by Bluetooth SIG
    RBLE_SVC_USB_ASSIGNED_ID            Vendor ID assigned by USB Implementer's
                                         Forum
};
```

Declaration of enumerated type for Name Space field of Characteristic Presentation Format descriptor

```
enum RBLE_SVC_PRESEN_NAMESPASE_enum {
    RBLE_SVC_NAMESPACE_SIG = 0x01,        Defined by Bluetooth SIG
};
```

Declaration of enumerated type for security level of Service

```
enum RBLE_SVC_SEC_LVL_enum {
    RBLE_SVC_SEC_NONE = 0x01,             No security
    RBLE_SVC_SEC_UNAUTH = 0x02,           Require unauthenticated pairing
    RBLE_SVC_SEC_AUTH = 0x04,             Require authenticated pairing
    RBLE_SVC_SEC_AUTZ = 0x08,             Require authorization
    RBLE_SVC_SEC_ENC = 0x10               Require encryption
};
```

Declaration of enumerated type for connection types

```
enum RBLE_PRF_CON_enum {
    RBLE_PRF_CON_DISCOVERY = 0x00,        Configuration connection performed
                                         when connecting for the first time
    RBLE_PRF_CON_NORMAL                    Normal connection performed when
                                         connecting for the second and
                                         subsequent times
};
```

Declaration of enumerated type for client configuration characteristic value

```
enum RBLE_PRF_CLIENT_CONFIG_enum {  
    RBLE_PRF_STOP_NTFFIND = 0x00,           Stop notification or indication of  
                                              characteristic value.  
    RBLE_PRF_START_NTF,                     Start notification of  
                                              characteristic value.  
    RBLE_PRF_START_IND                     Start indication of  
                                              characteristic value.  
};
```

Declaration of enumerated type for server configuration characteristic value

```
enum RBLE_PRF_SERVER_CONFIG_enum {  
    RBLE_PRF_STOP_BRD = 0x00,               Stop broadcast of characteristic value.  
    RBLE_PRF_START_BRD                     Start broadcast of characteristic value.  
};
```

## 2.2 Status Definitions

This section describes the status definitions used by the API of each profile.

- Declaration of enumerated type for rBLE status

```
enum RBLE_STATUS_enum {
    RBLE_OK                                = 0x00,    Normal operation
    RBLE_PRF_ERR_INVALID_PARAM            = 0x90,    Invalid parameter specified for
                                                    setting or acquiring a characteristic
                                                    value
    RBLE_PRF_ERR_INEXISTENT_HDL,          Invalid handle specified for setting
                                                    or acquiring a characteristic value
    RBLE_PRF_ERR_STOP_DISC_CHAR_MISSING,  The characteristic value is missing.
    RBLE_PRF_ERR_MULTIPLE_IAS,            Multiple IASs exist.
    RBLE_PRF_ERR_INCORRECT_PROP,          Incorrect property
    RBLE_PRF_ERR_MULTIPLE_CHAR,           Multiple characteristic values exist.
    RBLE_PRF_ERR_NOT_WRITABLE,            Writing is not permitted.
    RBLE_PRF_ERR_NOT_READABLE,           Reading is not permitted.
    RBLE_PRF_ERR_REQ_DISALLOWED,         Requesting is not permitted.
    RBLE_PRF_ERR_NTF_DISABLED,           Notification is disabled.
    RBLE_PRF_ERR_IND_DISABLED,           Indication is disabled.
    RBLE_PRF_ERR_ATT_NOT_SUPPORTED,       The characteristic value is not
                                                    supported.
};
```

Note: Statuses other than the above are described in *API Reference Manual: Basics*.

### 3. Blood Pressure Profile

This section describes the API of the Blood Pressure profile. The Blood Pressure profile is used to enable a device to obtain blood pressure measurement data from a blood pressure sensor.

#### 3.1 Definitions

This section describes the definitions used by the API of the Blood Pressure profile.

- Declaration of enumerated type for BLP event types

enum RBLE_BLP_EVENT_TYPE_enum {	
RBLE_BLP_EVENT_SENSOR_ENABLE_COMP = 0x01,	Sensor enable completion event (Parameter: sensor_enable)
RBLE_BLP_EVENT_SENSOR_DISABLE_COMP,	Sensor disable completion event (Parameter: sensor_disable)
RBLE_BLP_EVENT_SENSOR_ERROR_IND,	Sensor error indication event (Parameter: error_ind)
RBLE_BLP_EVENT_SENSOR_SEND_MEASUREMENTS_COMP,	Measured value send completion event (Parameter: send_measurements)
RBLE_BLP_EVENT_SENSOR_CFG_INDNTF_IND,	Sensor characteristic configuration change indication event (Parameter: blps_cfg_indntf_ind)
RBLE_BLP_EVENT_SENSOR_COMMAND_DISALLOWED_IND,	Command disallowed indication event (Parameter: cmd_disallowed_ind)
RBLE_BLP_EVENT_COLLECTOR_ENABLE_COMP = 0x81,	Collector enable completion event (Parameter: collector_enable)
RBLE_BLP_EVENT_COLLECTOR_DISABLE_COMP,	Collector disable completion event (Parameter: collector_disable)
RBLE_BLP_EVENT_COLLECTOR_ERROR_IND,	Collector error indication event (Parameter: error_ind)
RBLE_BLP_EVENT_COLLECTOR_MEASUREMENTS_IND,	Collector measured value indication event (Parameter: measurements_ind)
RBLE_BLP_EVENT_COLLECTOR_READ_CHAR_RESPONSE,	Characteristic value read request response event (Parameter: rd_char_resp)
RBLE_BLP_EVENT_COLLECTOR_WRITE_CHAR_RESPONSE,	Characteristic value write request response event (Parameter: wr_char_resp)

```
RBLE_BLP_EVENT_COLLECTOR_COMMAND_DISALLOWED_IND
```

```
Command disallowed indication
event
(Parameter: cmd_disallowed_ind)
```

```
};
```

- Declaration of data type for BLP event types

```
typedef uint8_t RBLE_BLP_EVENT_TYPE;
```

- Declaration of data type for BLP Sensor event callback function

```
typedef void ( *RBLE_BLPS_EVENT_HANDLER )( RBLE_BLPS_EVENT *event );
```

- Declaration of data type for BLP Collector event callback function

```
typedef void ( *RBLE_BLPC_EVENT_HANDLER )( RBLE_BLPC_EVENT *event );
```

- Declaration of enumerated type for blood pressure service/device information service characteristic codes

```
enum RBLE_BLPC_RD_CHAR_CODE_enum {
    RBLE_BLPC_RD_BLS_BM_CFG = 0x00,           Measurement result indication
    RBLE_BLPC_RD_BLS_IC_CFG,                  Intermediate cuff pressure
                                              information notification
    RBLE_BLPC_RD_BLS_BF,                      Blood pressure feature
    RBLE_BLPC_RD_DIS_MANUF,                   Blood pressure sensor system ID
    RBLE_BLPC_RD_DIS_MODEL,                   Blood pressure sensor model
                                              number
    RBLE_BLPC_RD_DIS_SERNB,                   Blood pressure sensor serial
                                              number
    RBLE_BLPC_RD_DIS_HWREV,                   Blood pressure sensor firmware
                                              revision
    RBLE_BLPC_RD_DIS_FWREV,                   Blood pressure sensor hardware
                                              revision
    RBLE_BLPC_RD_DIS_SWREV,                   Blood pressure sensor software
                                              revision
    RBLE_BLPC_RD_DIS_SYSID,                   Blood pressure sensor
                                              manufacturer name
    RBLE_BLPC_RD_DIS_IEEE,                   Blood pressure sensor IEEE
                                              certification information
};
```

- Declaration of enumerated type for blood pressure service characteristic value settings

```
enum RBLE_BLPC_WR_CHAR_CODE_enum {
    RBLE_BLPC_BLDPRS_MEAS_CODE = 0x01,       Measurement result indication
                                              setting
    RBLE_BLPC_INTERM_CUFPRS_CODE             Intermediate cuff pressure
                                              information notification setting
};
```

- Blood pressure service characteristic information structures

```
typedef struct RBLE_BLP_SENSOR_PARAM_t{
    uint16_t      bldprs_meas_ind_en;           Measurement result indication
                                                    configuration value
    uint16_t      interm_cufprs_ntf_en;         Intermediate cuff pressure information
                                                    notification configuration value
}RBLE_BLP_SENSOR_PARAM;
```

- Blood pressure measurement information structures

```
typedef struct RBLE_BLP_MEASUREMENTS_INFO_t {
    uint8_t        flag_stable_meas;           Measurement-in-progress flag
    uint8_t        flags;                     Data field flag
    int16_t        press_val_field1;           Systolic blood pressure value/
                                                    Cuff pressure
    int16_t        press_val_field2;           Diastolic blood pressure value/
                                                    Subfield 1
    int16_t        press_val_field3;           Average blood pressure/Subfield 2
    RBLE_DATE_TIME stamp;                     Time stamp
    int16_t        rate;                      Heart rate
    uint8_t        id;                       User ID
    uint8_t        reserved;                 Reserved
    uint16_t       meas_sts;                  Measurement status
}RBLE_BLP_MEASUREMENTS_INFO;
```



- Blood pressure service content structures

```
typedef struct RBLE_BLS_CONTENT_t{
    uint16_t      shdl;                Blood pressure service start handle
    uint16_t      ehdl;                Blood pressure service end handle
    uint16_t      bldprs_meas_char_hdl; Blood pressure measurement
                                      characteristic handle
    uint16_t      bldprs_meas_val_hdl;  Blood pressure measurement
                                      characteristic value handle
    uint16_t      bldprs_meas_cfg_hdl;  Blood pressure measurement client
                                      characteristic configuration
                                      descriptor handle
    uint8_t       bldprs_meas_prop;     Blood pressure measurement
                                      characteristic property
    uint8_t       reserved;             Reserved
    uint16_t      interm_cufprs_char_hdl; Intermediate cuff pressure
                                      characteristic handle
    uint16_t      interm_cufprs_val_hdl; Intermediate cuff pressure
                                      characteristic value handle
    uint16_t      interm_cufprs_cfg_hdl; Intermediate cuff pressure client
                                      characteristic configuration
                                      descriptor handle
    uint8_t       interm_cufprs_prop;   Intermediate cuff pressure
                                      characteristic property
    uint8_t       reserved2;           Reserved
    uint16_t      bldprs_feat_char_hdl; Blood pressure feature
                                      characteristic handle
    uint16_t      bldprs_feat_val_hdl;  Blood pressure feature
                                      characteristic value handle
    uint8_t       bldprs_feat_prop;     Blood pressure feature
                                      characteristic property
    uint8_t       reserved3;           Reserved
}RBLE_BLS_CONTENT;
```

- BLP Sensor event parameter structures

```
typedef struct RBLE_BLP_EVENT_t {
    RBLE_BLP_EVENT_TYPE    type;        BLP event type
    uint8_t                reserved;     Reserved
    union Event_Bls_Parameter_u {
        Generic event
        RBLE_STATUS        status;      Status

        Sensor enable completion event
        struct RBLE_BLP_Sensor_Enable_t{
            RBLE_STATUS    status;        Status
            uint8_t         reserved;     Reserved
            uint16_t        conhdl;       Connection handle
        }sensor_enable;
    };
}
```

**Sensor disable completion event**

```

struct RBLE_BLP_Sensor_Disable_t{
    uint16_t                conhdl;                Connection handle
    RBLE_BLP_SENSOR_PARAM   sensor_info;           Blood pressure service
                                                    information
}sensor_disable;

```

**Sensor error indication event**

```

struct RBLE_BLP_Sensor_Error_Ind_t{
    uint16_t                conhdl;                Connection handle
    RBLE_STATUS             status;                Status
}error_ind;

```

**Sensor measured value send completion event**

```

struct RBLE_BLP_Sensor_Send_Measurements_t{
    uint16_t                conhdl;                Connection handle
    RBLE_STATUS             status;                Status
}send_measurements;

```

**Sensor characteristic configuration change indication event**

```

struct RBLE_BLP_Sensor_Cfg_Indntf_Ind_t{
    uint16_t                conhdl;                Connection handle
    uint8_t                 char_code;             Status
    uint8_t                 reserved;              Reserved
    uint16_t                cfg_val;               Configuration
                                                    characteristic value
}blps_cfg_indntf_ind;

```

**Sensor command disallowed indication event**

```

struct RBLE_BLP_Sensor_Command_Disallowed_Ind_t{
    RBLE_STATUS             status;                Status
    uint8_t                 reserved;              Reserved
    uint16_t                opcode;               Opcode
}cmd_disallowed_ind;
} param;
} RBLE_BLPS_EVENT;

```

- BLP Collector event parameter structures

```
typedef struct RBLE_BLP_EVENT_t {
    RBLE_BLP_EVENT_TYPE      type;           BLP event type
    uint8_t                  reserved;       Reserved
    union Event_Blc_Parameter_u {
        Generic event
        RBLE_STATUS          status;        Status

        Collector enable completion event
        struct RBLE_BLP_Collector_Enable_t{
            RBLE_STATUS      status;        Status
            uint8_t          reserved;      Reserved
            uint16_t         conhdl;       Connection handle
            RBLE_BLS_CONTENT bls;          Blood pressure service
                                           content
            RBLE_DIS_CONTENT dis;          Device information service
                                           content
        }collector_enable;

        Collector disable completion event
        struct RBLE_BLP_Collector_Disable_t{
            RBLE_STATUS      status;        Status
            uint8_t          reserved;      Reserved
            uint16_t         conhdl;       Connection handle
        }collector_disable;

        Collector error indication event
        struct RBLE_BLP_Collector_Error_Ind_t{
            RBLE_STATUS      status;        Status
            uint8_t          reserved;      Reserved
            uint16_t         conhdl;       Connection handle
        }error_ind;

        Collector blood pressure measurement information indication event
        struct RBLE_BLP_Collector_Measurements_Ind_t{
            uint16_t         conhdl;       Connection handle
            RBLE_BLP_MEASUREMENTS_INFO measurements_info; Blood pressure
                                           measurement
                                           information
        }measurements_ind;

        Collector characteristic value read request response event
        struct RBLE_BLP_Collector_Read_Char_Response_t{
            uint16_t         conhdl;       Connection handle
            uint8_t          att_code;     Status
            RBLE_ATT_INFO_DATA data;       Acquired characteristic
                                           data
        }
    }
};
```

```
}rd_char_resp;
```

**Collector characteristic value write request response event**

```
struct RBLE_BLP_Collector_Write_Char_Response_t{
    uint16_t          conhdl;          Connection handle
    uint8_t           att_code;        Status
}wr_char_resp;
```

**Collector command disallowed indication event**

```
struct RBLE_BLP_Collector_Command_Disallowed_Ind_t{
    RBLE_STATUS       status;          Status
    uint8_t           reserved;        Reserved
    uint16_t          opcode;          Opcode
}cmd_disallowed_ind;
} param;
} RBLE_BLPC_EVENT;
```

## 3.2 Functions

The following table shows the API functions defined for the BLP of rBLE and the following sections describe the API functions in detail.

Table 3-1 API Functions Used by the BLP

RBLE_BLP_Sensor_Enable	Enables the Sensor role.
RBLE_BLP_Sensor_Disable	Disables the Sensor role.
RBLE_BLP_Sensor_Send_Measurements	Sends the measured value.
RBLE_BLP_Collector_Enable	Enables the Collector role.
RBLE_BLP_Collector_Disable	Disables the Collector role.
RBLE_BLP_Collector_Read_Char	Reads the characteristic value.
RBLE_BLP_Collector_Write_Char	Writes the characteristic value.

### 3.2.1 RBLE\_BLP\_Sensor\_Enable

**RBLE\_STATUS RBLE\_BLP\_Sensor\_Enable(uint16\_t conhdl, uint8\_t sec\_lv, uint8\_t con\_type, RBLE\_BLP\_SENSOR\_PARAM \*param, RBLE\_BLPS\_EVENT\_HANDLER call\_back)**

This function enables the BLP Sensor role.

If the measurement result indication and intermediate cuff pressure information notification setting has been specified from the Collector, set the indication/notification setting parameter to 0 to configure the connection. If this setting or information has been specified from the Sensor, perform a normal connection in accordance with the indication/notification setting parameter.

The result is reported by using the Sensor role enable completion event RBLE\_BLP\_EVENT\_SENSOR\_ENABLE\_COMP.

Parameters:

<i>conhdl</i>	Connection handle		
<i>sec_lv</i>	Security level		
<i>con_type</i>	RBLE_PRF_CON_DISCOVERY		Configuration connection
	RBLE_PRF_CON_NORMAL		Normal connection
<i>*param</i>	<i>bldprs_meas_ind_en</i>	RBLE_PRF_STOP_NTFFIND	Stop indication of the measurement result.
		RBLE_PRF_START_IND	Start indication of the measurement result.
	<i>interm_cufprs_ntf_en</i>	RBLE_PRF_STOP_NTFFIND	Stop notification of the intermediate cuff pressure information.
		RBLE_PRF_START_NTF	Start notification of the intermediate cuff pressure information.
<i>call_back</i>	Specify the callback function that reports the BLP event.		

Return:

<i>RBLE_OK</i>	Success
<i>RBLE_ERR</i>	Error occurred in BLP Sensor role enable processing
<i>RBLE_PARAM_ERR</i>	Invalid parameter
<i>RBLE_STATUS_ERROR</i>	Not executable because the rBLE mode is other than RBLE_MODE_ACTIVE.

### 3.2.2 RBLE\_BLP\_Sensor\_Disable

**RBLE\_STATUS RBLE\_BLP\_Sensor\_Disable(uint16\_t conhdl)**

This function disables the BLP Sensor role.

The result is reported by using the Sensor role disable completion event RBLE\_BLP\_EVENT\_SENSOR\_DISABLE\_COMP.

Parameters:

<i>conhdl</i>	Connection handle
---------------	-------------------

Return:

<i>RBLE_OK</i>	Success
<i>RBLE_STATUS_ERROR</i>	Not executable because the rBLE mode is other than RBLE_MODE_ACTIVE.

## 3.2.3 RBLE\_BLP\_Sensor\_Send\_Measurements

```
RBLE_STATUS RBLE_BLP_Sensor_Send_Measurements(uint16_t conhdl,
        RBLE_BLP_MEASUREMENTS_INFO *measurements_info)
```

This function sends the measured value data from the blood pressure sensor.

To send intermediate cuff pressure information from the blood pressure sensor, set `flag_stable_meas` to 0 and save the cuff pressure information in `press_val_field1` before executing this function. To send the measurement result after the blood pressure sensor has completed measuring the cuff pressure, set `flag_stable_meas` to 1 and save the cuff pressure information in `press_val_field1`, `press_val_field2`, and `press_val_field3` before executing this function.

The result is reported by using the Sensor role measured value send completion event `RBLE_BLP_EVENT_SENSOR_SEND_MEASUREMENTS_COMP`.

Parameters:

<i>conhdl</i>		Connection handle	
<i>*measurements_info</i>	<i>flag_stable_meas</i>		Flag indicating that measurement is in progress (0) or that measurement is not in progress (1)
	<i>flags</i>		Flag that defines whether there is a data field in the characteristic value or not
	<i>press_val_field1</i>		Measurement result: Systolic blood pressure value Measurement-in-progress: Cuff pressure
	<i>press_val_field2</i>		Measurement result: Diastolic blood pressure value Measurement-in-progress: Subfield 1 (normal value: 0x07FF)
	<i>press_val_field3</i>		Measurement result: Average blood pressure Measurement-in-progress: Subfield 2 (normal value: 0x07FF)
	<i>stamp</i>	<i>year</i>	Year
		<i>month</i>	Month
		<i>day</i>	Day
		<i>hour</i>	Hour
		<i>min</i>	Minute
		<i>sec</i>	Second
	<i>rate</i>		Heart rate
	<i>id</i>		User ID
	<i>meas_sts</i>		Measurement status

Return:

<i>RBLE_OK</i>	Success
<i>RBLE_STATUS_ERROR</i>	Not executable because the rBLE mode is other than <code>RBLE_MODE_ACTIVE</code> .

## 3.2.4 RBLE\_BLP\_Collector\_Enable

```
RBLE_STATUS RBLE_BLP_Collector_Enable(uint16_t conhdl, uint8_t con_type, RBLE_BLS_CONTENT *bls,
RBLE_DIS_CONTENT *dis, RBLE_BLPC_EVENT_HANDLER call_back)
```

This function enables the BLP Collector role and starts access to the services exposed by the BLP Sensor. The result is reported by using the Collector role enable completion event

RBLE\_BLP\_EVENT\_COLLECTOR\_ENABLE\_COMP.

When starting access to the service exposed by a BLP Sensor to be connected for the first time, set 0 to the parameter of the service to configure the connection and to discover the service for the Sensor. If the handle information about the discovered service is saved and is used when the Sensor is connected for a second or subsequent time, detecting the service is skipped, which enables a high-speed access to the service.

While the Collector role is enabled, the service exposed by only one BLP Sensor is accessible. To connect to more than one Sensor at the same time and access the services exposed by each Sensor, repeat enable/disable of the Collector role in order to switch access to them. At that time, perform normal connection by using the connection handle (which was obtained when connecting to each Sensor) and the handle information (which was saved when starting access to the service for the first time) as parameters.

Parameters:

<i>conhdl</i>	Connection handle	
<i>con_type</i>	RBLE_PRF_CON_DISCOVERY	Configuration connection performed when connecting for the first time
	RBLE_PRF_CON_NORMAL	Normal connection performed when connecting for the second and subsequent times
<i>*bls</i>	<i>shdl</i>	Blood pressure service start handle
	<i>ehdl</i>	Blood pressure service end handle
	<i>bldprs_meas_char_hdl</i>	Blood pressure measurement characteristic handle
	<i>bldprs_meas_val_hdl</i>	Blood pressure measurement characteristic value handle
	<i>bldprs_meas_cfg_hdl</i>	Blood pressure measurement client characteristic configuration descriptor handle
	<i>bldprs_meas_prop</i>	Blood pressure measurement characteristic property
	<i>interm_cufprs_char_hdl</i>	Intermediate cuff pressure characteristic handle
	<i>interm_cufprs_val_hdl</i>	Intermediate cuff pressure characteristic value handle
	<i>interm_cufprs_cfg_hdl</i>	Intermediate cuff pressure client characteristic configuration descriptor handle
	<i>interm_cufprs_prop</i>	Intermediate cuff pressure characteristic property
	<i>bldprs_feat_char_hdl</i>	Blood pressure feature characteristic handle
	<i>bldprs_feat_val_hdl</i>	Blood pressure feature characteristic value handle
	<i>bldprs_feat_prop</i>	Blood pressure feature characteristic property
<i>*dis</i>	<i>shdl</i>	Device information service start handle
	<i>ehdl</i>	Device information service end handle
	<i>sys_id_char_hdl</i>	System ID characteristic handle
	<i>sys_id_val_hdl</i>	System ID characteristic value handle
	<i>sys_id_prop</i>	System ID characteristic property
	<i>model_nb_char_hdl</i>	Model number characteristic handle
	<i>model_nb_val_hdl</i>	Model number characteristic value handle
	<i>model_nb_prop</i>	Model number characteristic property



	<i>serial_nb_char_hdl</i>	Serial number characteristic handle
	<i>serial_nb_val_hdl</i>	Serial number characteristic value handle
	<i>serial_nb_prop</i>	Serial number characteristic property
	<i>fw_rev_char_hdl</i>	Firmware revision characteristic handle
	<i>fw_rev_val_hdl</i>	Firmware revision characteristic value handle
	<i>fw_rev_prop</i>	Firmware revision characteristic property
	<i>hw_rev_char_hdl</i>	Hardware revision characteristic handle
	<i>hw_rev_val_hdl</i>	Hardware revision characteristic value handle
	<i>hw_rev_prop</i>	Hardware revision characteristic property
	<i>sw_rev_char_hdl</i>	Software revision characteristic handle
	<i>sw_rev_val_hdl</i>	Software revision characteristic value handle
	<i>sw_rev_prop</i>	Software revision characteristic property
	<i>manuf_name_char_hdl</i>	Manufacturer name characteristic handle
	<i>manuf_name_val_hdl</i>	Manufacturer name characteristic value handle
	<i>manuf_name_prop</i>	Manufacturer name characteristic property
	<i>ieee_certif_char_hdl</i>	IEEE certification characteristic handle
	<i>ieee_certif_val_hdl</i>	IEEE certification characteristic value handle
	<i>ieee_certif_prop</i>	IEEE certification characteristic property
	<i>call_back</i>	Specify the callback function that reports the BLP event.
Return:		
	<i>RBLE_OK</i>	Success
	<i>RBLE_ERR</i>	Error occurred in BLP Collector role enable processing
	<i>RBLE_PARAM_ERR</i>	Invalid parameter
	<i>RBLE_STATUS_ERROR</i>	Not executable because the rBLE mode is other than RBLE_MODE_ACTIVE.

### 3.2.5 RBLE\_BLP\_Collector\_Disable

<b>RBLE_STATUS RBLE_BLP_Collector_Disable(uint16_t conhdl)</b>		
This function disables the BLP Collector role and terminates access to the service exposed by the BLP Sensor. The result is reported by using the Collector role disable completion event RBLE_BLP_EVENT_COLLECTOR_DISABLE_COMP.		
Parameters:		
	<i>conhdl</i>	Connection handle
Return:		
	<i>RBLE_OK</i>	Success
	<i>RBLE_STATUS_ERROR</i>	Not executable because the rBLE mode is other than RBLE_MODE_ACTIVE.

## 3.2.6 RBLE\_BLP\_Collector\_Read\_Char

**RBLE\_STATUS RBLE\_BLP\_Collector\_Read\_Char(uint16\_t conhdl, uint8\_t char\_code)**

This function reads the characteristic value of the blood pressure service and the device information service.  
The result is reported by using the characteristic value read request response event  
RBLE\_BLP\_EVENT\_COLLECTOR\_READ\_CHAR\_RESPONSE.

Parameters:

<i>conhdl</i>	Connection handle	
<i>char_code</i>	RBLE_BLPC_RD_BLS_BM_CFG	Measurement result indication setting
	RBLE_BLPC_RD_BLS_IC_CFG	Intermediate cuff pressure information notification setting
	RBLE_BLPC_RD_BLS_BF	Blood pressure feature setting
	RBLE_BLPC_RD_DIS_SYSID	Blood pressure sensor system ID
	RBLE_BLPC_RD_DIS_MODEL	Blood pressure sensor model number
	RBLE_BLPC_RD_DIS_SERNB	Blood pressure sensor serial number
	RBLE_BLPC_RD_DIS_FWREV	Blood pressure sensor firmware revision
	RBLE_BLPC_RD_DIS_HWREV	Blood pressure sensor hardware revision
	RBLE_BLPC_RD_DIS_SWREV	Blood pressure sensor software revision
	RBLE_BLPC_RD_DIS_MANUF	Blood pressure sensor manufacturer name
	RBLE_BLPC_RD_DIS_IEEE	Blood pressure sensor IEEE certification information

Return:

<i>RBLE_OK</i>	Success
<i>RBLE_STATUS_ERROR</i>	Not executable because the rBLE mode is other than RBLE_MODE_ACTIVE.

## 3.2.7 RBLE\_BLP\_Collector\_Write\_Char

**RBLE\_STATUS RBLE\_BLP\_Collector\_Write\_Char(uint16\_t conhdl, uint8\_t char\_code, uint16\_t cfg\_val)**

This function writes the characteristic configuration descriptor of the blood pressure service and the device information service.

The result is reported by using the characteristic value write request response event  
RBLE\_BLP\_EVENT\_COLLECTOR\_WRITE\_CHAR\_RESPONSE.

Parameters:

<i>conhdl</i>	Connection handle	
<i>char_code</i>	RBLE_BLPC_BLDPRS_MEAS_CODE	Measurement result indication setting
	RBLE_BLPC_INTERM_CUFPRS_CODE	Intermediate cuff pressure information notification setting
<i>cfg_val</i>	RBLE_PRF_STOP_NTFFIND	Stop notification or indication.
	RBLE_PRF_START_NTF	Start notification.
	RBLE_PRF_START_IND	Start indication.

Return:

<i>RBLE_OK</i>	Success
<i>RBLE_STATUS_ERROR</i>	Not executable because the rBLE mode is other than RBLE_MODE_ACTIVE.

### 3.3 Events

The following table shows the events defined for the BLP of rBLE and the following sections describe the events in detail.

Table 3-2 Events Defined for the BLP

RBLE_BLP_EVENT_SENSOR_ENABLE_COMP	Sensor role enable completion event
RBLE_BLP_EVENT_SENSOR_DISABLE_COMP	Sensor role disable completion event
RBLE_BLP_EVENT_SENSOR_ERROR_IND	Sensor role error indication event
RBLE_BLP_EVENT_SENSOR_SEND_MEASUREMENTS_COMP	Measured value send completion event
RBLE_BLP_EVENT_SENSOR_CFG_INDNTF_IND	Characteristic value indication event
RBLE_BLP_EVENT_SENSOR_COMMAND_DISALLOWED_IND	Sensor role command disallowed indication event
RBLE_BLP_EVENT_COLLECTOR_ENABLE_COMP	Collector role enable completion event
RBLE_BLP_EVENT_COLLECTOR_DISABLE_COMP	Collector role disable completion event
RBLE_BLP_EVENT_COLLECTOR_ERROR_IND	Collector role error indication event
RBLE_BLP_EVENT_COLLECTOR_MEASUREMENTS_IND	Measured value indication event
RBLE_BLP_EVENT_COLLECTOR_READ_CHAR_RESPONSE	Characteristic value read request response event
RBLE_BLP_EVENT_COLLECTOR_WRITE_CHAR_RESPONSE	Characteristic value write request response event
RBLE_BLP_EVENT_COLLECTOR_COMMAND_DISALLOWED_IND	Collector role command disallowed indication event

## 3.3.1 RBLE\_BLP\_EVENT\_SENSOR\_ENABLE\_COMP

RBLE_BLP_EVENT_SENSOR_ENABLE_COMP	
This event reports the result of enabling the BLP Sensor role (RBLE_BLP_Sensor_Role_Enable).	
Parameters:	
<i>status</i>	Result of enabling the Sensor role (See 2.2 and <i>Bluetooth Low Energy Protocol Stack API Reference Manual: Basics, 3.2, Declaration of enumerated type for rBLE status.</i> )
<i>conhdl</i>	Connection handle

## 3.3.2 RBLE\_BLP\_EVENT\_SENSOR\_DISABLE\_COMP

RBLE_BLP_EVENT_SENSOR_DISABLE_COMP			
This event reports the result of disabling the BLP Sensor role (RBLE_BLP_Sensor_Role_Disable).			
Parameters:			
<i>conhdl</i>	Connection handle		
<i>sensor_info</i>	<i>bldprs_meas_ind_en</i>	RBLE_PRF_STOP_NTFFIND	Stop notification/indication of the measurement result.
		RBLE_PRF_START_IND	Start indication of the measurement result.
	<i>interm_cufprs_ntf_en</i>	RBLE_PRF_STOP_NTFFIND	Stop notification/indication of the intermediate cuff pressure information.
		RBLE_PRF_START_NTF	Start notification of the intermediate cuff pressure information.

## 3.3.3 RBLE\_BLP\_EVENT\_SENSOR\_ERROR\_IND

RBLE_BLP_EVENT_SENSOR_ERROR_IND	
This event indicates an error code unique to the BLP Sensor role.	
Parameters:	
<i>conhdl</i>	Connection handle
<i>status</i>	Error code (See 2.2 and <i>Bluetooth Low Energy Protocol Stack API Reference Manual: Basics, 3.2, Declaration of enumerated type for rBLE status.</i> )

## 3.3.4 RBLE\_BLP\_EVENT\_SENSOR\_SEND\_MEASUREMENTS\_COMP

RBLE_BLP_EVENT_SENSOR_SEND_MEASUREMENTS_COMP	
This event reports completion of sending the measured value (RBLE_BLP_Sensor_Send_Measurements).	
Parameters:	
<i>conhdl</i>	Connection handle
<i>status</i>	Measured value send completion result (See 2.2 and <i>Bluetooth Low Energy Protocol Stack API Reference Manual: Basics, 3.2, Declaration of enumerated type for rBLE status.</i> )

## 3.3.5 RBLE\_BLP\_EVENT\_SENSOR\_CFG\_INDNTF\_IND

RBLE_BLP_EVENT_SENSOR_CFG_INDNTF_IND			
This event reports that the value of the client characteristic configuration descriptor of the blood pressure service has been written.			
Parameters:			
<i>conhdl</i>	Connection handle		
<i>char_code</i>	RBLE_BLPC_BLDPRS_MEAS_CODE	Measurement result indication setting	
	RBLE_BLPC_INTERM_CUFPRS_CODE	Intermediate cuff pressure information notification setting	
<i>cfg_val</i>	RBLE_PRF_STOP_NTFFIND	Stop notification or indication.	
	RBLE_PRF_START_NTF	Start notification.	
	RBLE_PRF_START_IND	Start indication.	

## 3.3.6 RBLE\_BLP\_EVENT\_SENSOR\_COMMAND\_DISALLOWED\_IND

RBLE_BLP_EVENT_SENSOR_COMMAND_DISALLOWED_IND			
This event indicates the error that occurs when a command executed by the Sensor role cannot be accepted.			
Parameters:			
<i>status</i>	Result of command execution (See 2.2 and <i>Bluetooth Low Energy Protocol Stack API Reference Manual: Basics, 3.2, Declaration of enumerated type for rBLE status.</i> )		
<i>opcode</i>	RBLE_CMD_BLP_SENSOR_ENABLE	Sensor enable command	
	RBLE_CMD_BLP_SENSOR_DISABLE	Sensor disable command	
	RBLE_CMD_BLP_SENSOR_SEND_MEASUREMENTS	Measured data send command	

## 3.3.7 RBLE\_BLP\_EVENT\_COLLECTOR\_ENABLE\_COMP

## RBLE\_BLP\_EVENT\_COLLECTOR\_ENABLE\_COMP

This event reports the result of enabling the BLP Collector role (RBLE\_BLP\_Collector\_Role\_Enable). Save the obtained handle information about the discovered service, to enable a high-speed access to the service without service detection when restarting access to the service.

Parameters:

<i>status</i>	Result of enabling the Collector role (See 2.2 and <i>Bluetooth Low Energy Protocol Stack API Reference Manual: Basics, 3.2, Declaration of enumerated type for rBLE status.</i> )	
	Connection handle	
<i>bls</i>	<i>shdl</i>	Blood pressure service start handle
	<i>ehdl</i>	Blood pressure service end handle
	<i>bldprs_meas_char_hdl</i>	Blood pressure measurement characteristic handle
	<i>bldprs_meas_val_hdl</i>	Blood pressure measurement characteristic value handle
	<i>bldprs_meas_cfg_hdl</i>	Blood pressure measurement client characteristic configuration descriptor handle
	<i>bldprs_meas_prop</i>	Blood pressure measurement characteristic property
	<i>interm_cufprs_char_hdl</i>	Intermediate cuff pressure characteristic handle
	<i>interm_cufprs_val_hdl</i>	Intermediate cuff pressure characteristic value handle
	<i>interm_cufprs_cfg_hdl</i>	Intermediate cuff pressure client characteristic configuration descriptor handle
	<i>interm_cufprs_prop</i>	Intermediate cuff pressure characteristic property
	<i>bldprs_feat_char_hdl</i>	Blood pressure feature characteristic handle
	<i>bldprs_feat_val_hdl</i>	Blood pressure feature characteristic value handle
	<i>bldprs_feat_prop</i>	Blood pressure feature characteristic property
<i>dis</i>	<i>shdl</i>	Start handle
	<i>ehdl</i>	End handle
	<i>sys_id_char_hdl</i>	System ID characteristic handle
	<i>sys_id_val_hdl</i>	System ID characteristic value handle
	<i>sys_id_prop</i>	System ID characteristic property
	<i>model_nb_char_hdl</i>	Model number characteristic handle
	<i>model_nb_val_hdl</i>	Model number characteristic value handle
	<i>model_nb_prop</i>	Model number characteristic property
	<i>serial_nb_char_hdl</i>	Serial number characteristic handle
	<i>serial_nb_val_hdl</i>	Serial number characteristic value handle
	<i>serial_nb_prop</i>	Serial number characteristic property
	<i>fw_rev_char_hdl</i>	Firmware revision characteristic handle
	<i>fw_rev_val_hdl</i>	Firmware revision characteristic value handle
	<i>fw_rev_prop</i>	Firmware revision characteristic property
	<i>hw_rev_char_hdl</i>	Hardware revision characteristic handle
	<i>hw_rev_val_hdl</i>	Hardware revision characteristic value handle
	<i>hw_rev_prop</i>	Hardware revision characteristic property
	<i>sw_rev_char_hdl</i>	Software revision characteristic handle
	<i>sw_rev_val_hdl</i>	Software revision characteristic value handle
	<i>sw_rev_prop</i>	Software revision characteristic property
	<i>manuf_name_char_hdl</i>	Manufacturer name characteristic handle

		<i>manuf_name_val_hdl</i>	Manufacturer name characteristic value handle
		<i>manuf_name_prop</i>	Manufacturer name characteristic property
		<i>ieee_certif_char_hdl</i>	IEEE certification characteristic handle
		<i>ieee_certif_val_hdl</i>	IEEE certification characteristic value handle
		<i>ieee_certif_prop</i>	IEEE certification characteristic property

### 3.3.8 RBLE\_BLP\_EVENT\_COLLECTOR\_DISABLE\_COMP

RBLE_BLP_EVENT_COLLECTOR_DISABLE_COMP		
This event reports the result of disabling the BLP Collector role (RBLE_BLP_Collector_Role_Disable).		
Parameters:		
	<i>status</i>	Result of disabling the Collector role (See 2.2 and <i>Bluetooth Low Energy Protocol Stack API Reference Manual: Basics, 3.2, Declaration of enumerated type for rBLE status.</i> )
	<i>conhdl</i>	Connection handle

### 3.3.9 RBLE\_BLP\_EVENT\_COLLECTOR\_ERROR\_IND

RBLE_BLP_EVENT_COLLECTOR_ERROR_IND		
This event indicates an error code unique to the BLP Collector role.		
Parameters:		
	<i>status</i>	Error code (See 2.2 and <i>Bluetooth Low Energy Protocol Stack API Reference Manual: Basics, 3.2, Declaration of enumerated type for rBLE status.</i> )
	<i>conhdl</i>	Connection handle

## 3.3.10 RBLE\_BLP\_EVENT\_COLLECTOR\_MEASUREMENTS\_IND

## RBLE\_BLP\_EVENT\_COLLECTOR\_MEASUREMENTS\_IND

This event indicates the measured value sent from the Sensor.

When sending intermediate cuff pressure information, the cuff pressure information is indicated with `flag_stable_meas` set to 0 and then saved in `press_val_field1`. When sending the cuff pressure information after measurement is complete, the cuff pressure information is indicated with `flag_stable_meas` set to 1 and then saved in `press_val_field1`, `press_val_field2`, and `press_val_field3`.

Parameters:

<i>conhdl</i>	Connection handle	
	<i>flag_stable_meas</i>	Flag indicating that measurement is in progress (0) or that measurement is not in progress (1)
	<i>flags</i>	Flag that defines whether there is a data field in the characteristic value or not
	<i>press_val_field1</i>	Measurement result: Systolic blood pressure value Measurement-in-progress: Cuff pressure
	<i>press_val_field2</i>	Measurement result: Diastolic blood pressure value Measurement-in-progress: Subfield 1 (normal value: 0x07FF)
	<i>press_val_field3</i>	Measurement result: Average blood pressure Measurement-in-progress: Subfield 2 (normal value: 0x07FF)
	<i>stamp</i>	<i>year</i> Year
		<i>month</i> Month
		<i>day</i> Day
		<i>hour</i> Hour
		<i>min</i> Minute
		<i>sec</i> Second
	<i>rate</i>	Heart rate
	<i>id</i>	User ID
	<i>meas_sts</i>	Measurement status

## 3.3.11 RBLE\_BLP\_EVENT\_COLLECTOR\_READ\_CHAR\_RESPONSE

## RBLE\_BLP\_EVENT\_COLLECTOR\_READ\_CHAR\_RESPONSE

This event indicates the response to the characteristic value acquisition request (RBLE\_BLP\_Collector\_Read\_Char).

Read out the acquired data in accordance with the contents of the request.

Parameters:

<i>conhdl</i>	Connection handle	
<i>att_code</i>	0x00	Characteristic value successfully acquired
	Other than 0x00	Error occurred when acquiring characteristic value
<i>data</i>	<i>each_len</i>	Length of each result
	<i>len</i>	Data length
	<i>data[RBLE_ATT_MAX_VALUE]</i>	Read characteristic data



## 3.3.12 RBLE\_BLP\_EVENT\_COLLECTOR\_WRITE\_CHAR\_RESPONSE

RBLE_BLP_EVENT_COLLECTOR_WRITE_CHAR_RESPONSE			
This event indicates the response to the characteristic value write request (RBLE_BLP_Collector_Write_Char).			
Parameters:			
<i>conhdl</i>	Connection handle		
<i>att_code</i>	0x00	Characteristic value successfully written	
	Other than 0x00	Error occurred when writing characteristic value	

## 3.3.13 RBLE\_BLP\_EVENT\_COLLECTOR\_COMMAND\_DISALLOWED\_IND

RBLE_BLP_EVENT_COLLECTOR_COMMAND_DISALLOWED_IND			
This event indicates the error that occurs when a command executed by the Collector role cannot be accepted.			
Parameters:			
<i>status</i>	Result of command execution (See 2.2 and <i>Bluetooth Low Energy Protocol Stack API Reference Manual: Basics, 3.2, Declaration of enumerated type for rBLE status.</i> )		
<i>opcode</i>	RBLE_CMD_BLP_COLLECTOR_ENABLE	Collector enable command	
	RBLE_CMD_BLP_COLLECTOR_DISABLE	Collector disable command	
	RBLE_CMD_BLP_COLLECTOR_READ_CHAR	Characteristic read command	
	RBLE_CMD_BLP_COLLECTOR_WRITE_CHAR	Characteristic write command	

## 3.4 Message Sequence Chart




Figure 3-1 Example of Use Case In Which BLP Is Implemented by Using rBLE API

## 4. Notes

## Appendix A How to Read Definition Tables

This section shows how to read the tables that describes the rBLE API functions and events shown in this document.

### A.1 How to Read Function Definition Tables

The following contents are included in the function definition tables:

The diagram illustrates the structure of a function definition table. It consists of several sections: a light green header row for the function prototype, a row for the function operation and event, a 'Parameters' section with a table of parameter details, and a 'Return' section with a table of return values. Callouts provide detailed explanations for each part.

**Callout 1:** The Parameters area describes the parameters specified for the function. The italicized character strings on the left are the parameters of the function. The meaning of each parameter is described on the far right following the variables.

**Callout 2:** The italicized character string(s) next to each parameter indicate the member(s) of the parameter (structure).

**Callout 3:** The values that can be specified for the parameter might be described between the parameter name and its description.

**Table Structure:**

The function definition is shown at the top of the table in the row with the light green background. This area shows the function prototype.				
The operation of the function and the event reported after executing the function are described in this area.				
Parameters:				
<i>Parameter 1</i>	Description of parameter 1			
<i>Parameter 2</i>	<i>Member 1</i>	Value 1 that can be specified for member 1	Description of value 1 that can be specified for member 1	
		Value 1 that can be specified for member 2	Description of value 1 that can be specified for member 2	
	<i>Member 2</i>	Description of member 2		
Return:				
<i>Value 1 that might be returned</i>		Description of value 1 that might be returned		
<i>Value 2 that might be returned</i>		Description of value 2 that might be returned		

**Callout 4:** The Return area describes the values returned for the function. The leftmost row shows the value that might be returned, and the next row describes the return value.

## A.2 How to Read Event Definition Tables

The following contents are included in the event definition tables:

The diagram illustrates the structure of an event definition table. It consists of several rows and columns. The top row has an orange background and contains the event definition. The second row describes the information reported by the event. The third row is labeled 'Parameters:' and contains a table of parameters. The parameter table has three main rows: 'Parameter 1', 'Parameter 2', and 'Parameter 3'. 'Parameter 1' has a single description. 'Parameter 2' has three members, each with a description. 'Parameter 3' has two values, each with a description. Callouts explain the italicized character strings and the values specified for the parameters.

The event definition is shown at the top of the table in the row with the orange background. This area shows the event type.		
The information reported by the event is described in this area.		
Parameters:		
<i>Parameter 1</i>	Description of parameter 1	
<i>Parameter 2</i>	<i>Member 1</i>	Description of member 1
	<i>Member 2</i>	Description of member 2
	<i>Member 3</i>	Description of member 3
<i>Parameter 3</i>	Value 1 that can be specified for parameter 3	Description of value 1 that can be specified for parameter 3
	Value 2 that can be specified for parameter 3	Description of value 2 that can be specified for parameter 3

The Parameters area describes the parameters specified for the event. The italicized character strings on the left show the parameters of the event parameter structure. The meaning of each parameter is described on the far right.

The italicized character string(s) next to each parameter indicate the member(s) of the parameter (structure).

The values that can be specified for the parameter might be shown between the parameter name and its description.

## Appendix B Referenced Documents

1. Bluetooth Core Specification v4.0, Bluetooth SIG
2. Find Me Profile Specification v1.0, Bluetooth SIG
3. Immediate Alert Service Specification v1.0, Bluetooth SIG
4. Proximity Profile Specification v1.0, Bluetooth SIG
5. Link Loss Service Specification v1.0, Bluetooth SIG
6. Tx Power Service Specification v1.0, Bluetooth SIG
7. Health Thermometer Profile Specification v1.0, Bluetooth SIG
8. Health Thermometer Service Specification v1.0, Bluetooth SIG
9. Device Information Service Specification v1.1, Bluetooth SIG
10. Blood Pressure Profile Specification v1.0, Bluetooth SIG
11. Blood Pressure Service Specification v1.0, Bluetooth SIG
12. HID over GATT Profile Specification v1.0, Bluetooth SIG
13. HID Service Specification v1.0, Bluetooth SIG
14. Battery Service Specification v1.0, Bluetooth SIG
15. Scan Parameters Profile Specification v1.0, Bluetooth SIG
16. Scan Parameters Service Specification v1.0, Bluetooth SIG
17. Bluetooth SIG Assigned Numbers <https://www.bluetooth.org/Technical/AssignedNumbers/home.htm>
18. Services & Characteristics UUID <http://developer.bluetooth.org/gatt/Pages/default.aspx>
19. Personal Health Devices Transcoding White Paper v1.2, Bluetooth SIG

## Appendix C Terminology

Term	Description
Service	A service is provided from a GATT server to a GATT client. The GATT server exposes some characteristics as the interface. The service prescribes how to access the exposed characteristics.
Profile	A profile enables implementation of a use case by using one or more services. The services used are defined in the specifications of each profile.
Characteristic	A characteristic is a value used to identify services. The characteristics to be exposed and their formats are defined by each service.
Role	Each device takes the role prescribed by the profile or service in order to implement the specified use case.
Client Characteristic Configuration Descriptor	A descriptor is used to control notifications or indications of characteristic values that include the client characteristic configuration descriptor sent from the GATT server.
Connection Handle	The handle determined by the controller stack and is used to identify connection with a remote device. The valid handle range is between 0x0000 and 0x0EFF.

REVISION HISTORY	Bluetooth Low Energy Protocol Stack API Reference Manual: BLP
------------------	---

Rev.	Date	Description	
		Page	Summary
1.00	Feb 15, 2013	---	First Edition issued
1.01	Mar 27, 2013	---	The description about the high-speed access to the service for a second or subsequent time is added.
1.02	Jun 28, 2013	---	Bookmark is added.
1.03	Sep 19, 2014	2	The common definitions of profile are added.
		5	Definitions of client configuration characteristic value and connection type are deleted.
		---	Parameter description is changed to use the common definitions of profile.
1.04	Apr 17, 2015	2	The service definitions are updated.

a



---

Bluetooth Low Energy Protocol Stack  
API Reference Manual: BLP

Publication Date: Rev.1.04 Apr 17, 2015

Published by: Renesas Electronics Corporation

---



## SALES OFFICES

## Renesas Electronics Corporation

<http://www.renesas.com>

Refer to "<http://www.renesas.com/>" for the latest and detailed information.

### **Renesas Electronics America Inc.**

2801 Scott Boulevard Santa Clara, CA 95050-2549, U.S.A.  
Tel: +1-408-588-6000, Fax: +1-408-588-6130

### **Renesas Electronics Canada Limited**

9251 Yonge Street, Suite 8309 Richmond Hill, Ontario Canada L4C 9T3  
Tel: +1-905-237-2004

### **Renesas Electronics Europe Limited**

Dukes Meadow, Millboard Road, Bourne End, Buckinghamshire, SL8 5FH, U.K.  
Tel: +44-1628-585-100, Fax: +44-1628-585-900

### **Renesas Electronics Europe GmbH**

Arcadiastrasse 10, 40472 Düsseldorf, Germany  
Tel: +49-211-6503-0, Fax: +49-211-6503-1327

### **Renesas Electronics (China) Co., Ltd.**

Room 1709, Quantum Plaza, No.27 ZhiChunLu Haidian District, Beijing 100191, P.R.China  
Tel: +86-10-8235-1155, Fax: +86-10-8235-7679

### **Renesas Electronics (Shanghai) Co., Ltd.**

Unit 301, Tower A, Central Towers, 555 Langao Road, Putuo District, Shanghai, P. R. China 200333  
Tel: +86-21-2226-0888, Fax: +86-21-2226-0999

### **Renesas Electronics Hong Kong Limited**

Unit 1601-1611, 16/F., Tower 2, Grand Century Place, 193 Prince Edward Road West, Mongkok, Kowloon, Hong Kong  
Tel: +852-2265-6688, Fax: +852 2886-9022

### **Renesas Electronics Taiwan Co., Ltd.**

13F, No. 363, Fu Shing North Road, Taipei 10543, Taiwan  
Tel: +886-2-8175-9600, Fax: +886 2-8175-9670

### **Renesas Electronics Singapore Pte. Ltd.**

80 Bendemeer Road, Unit #06-02 Hyflux Innovation Centre, Singapore 339949  
Tel: +65-6213-0200, Fax: +65-6213-0300

### **Renesas Electronics Malaysia Sdn.Bhd.**

Unit 1207, Block B, Menara Amcorp, Amcorp Trade Centre, No. 18, Jln Persiaran Barat, 46050 Petaling Jaya, Selangor Darul Ehsan, Malaysia  
Tel: +60-3-7955-9390, Fax: +60-3-7955-9510

### **Renesas Electronics India Pvt. Ltd.**

No.777C, 100 Feet Road, HALII Stage, Indiranagar, Bangalore, India  
Tel: +91-80-67208700, Fax: +91-80-67208777

### **Renesas Electronics Korea Co., Ltd.**

12F., 234 Teheran-ro, Gangnam-Gu, Seoul, 135-080, Korea  
Tel: +82-2-558-3737, Fax: +82-2-558-5141

# Bluetooth Low Energy Protocol Stack



Renesas Electronics Corporation

R01UW0092EJ0104